Backstop in the right position during injection.

## en INSTRUCTIONS FOR USE FOR BELOTERO® BALANCE LIDOCAINE

### Description

BELOTERO Balance Lidocaine is a sterile, non-pyrogenic, viscoelastic, colourless, transparent gel of cross-linked sodium hyaluronate of non-animal origin, in a physiological phosphate buffer. BELOTERO Balance Lidocaine contains 0.3% of lidocaine hydrochloride.

### Presentation

BELOTERO Balance Lidocaine is presented in a 1 ml sterile glass syringe for single use, sterilized in autoclave by moist heat. Each box contains two sterile CE-marked needles of 27G½ (0.4x13mm) and 30G½ (0.3x13 mm) and two traceability labels per syringe as well as instructions for use.

### Composition

Cross-linked sodium hyaluronate: 22.5mg/ml  
Lidocaine hydrochloride: 3.0 mg/ml  
Phosphate buffer solution q.s.p.: 1.0 ml

The gel and all packaging items in direct contact with the gel are latex-free, bisphenol-A-free and phthalates-free. The volume of the gel in each syringe is stated on the external box.

### Indications

BELOTERO Balance Lidocaine is indicated for use to increase the volume of skin tissue which enables:

- Correction of facial depressions, either due to injury or age-related (scarring, moderate or deep wrinkles, furrows),
- Remodeling the contours of the face by injections into facial areas such as the cheeks, chin, etc.
- Increasing the volume of the lips (the border of the red part of the lips, or increasing the roundness)

The presence of lidocaine reduces local pain associated to the injection of the gel.

### Posology and administration method

This device is designed to be injected into the middle dermis by a legally approved practitioner. For successful treatment it is essential that the practitioner has received prior training in the injection technique for increasing the volume of skin tissues. The treatment must be carried out under the appropriate aseptic conditions. This product has been designed for use in physician's office.

BELOTERO Balance Lidocaine must be injected into a healthy, non-inflamed, disinfected skin. To ensure optimal use of BELOTERO Balance Lidocaine, it is recommended to assemble the needle according to the diagrams below. Start the injection pointing the needle bevel downwards.

The practitioner should note that the use of 30G½ needles will require greater force when injecting the gel than when using a 27G½.

BELOTERO Balance Lidocaine may be injected into the middle dermis using serial punctual injections, the linear retracing technique, anterotracing (push ahead) technique, cross-hatching technique (blanket), fan distribution technique, multi-puncture and micro-papular injection techniques. BELOTERO Balance Lidocaine may be used in combination with other Merz fillers respecting the correct depth of injection and indication for each product.

If the needle becomes obstructed and the injection pressure becomes too high, stop the injection and change the needle. Inject BELOTERO Balance Lidocaine slowly and apply the least amount of pressure necessary. The quantity to inject depends on the correction to be achieved. Over-correction is not necessary. After the injection, the practitioner can massage gently in order to distribute the product uniformly.

### Contra-indications

BELOTERO Balance Lidocaine is contra-indicated in case of:

- known hypersensitivity to one of the product's components, especially sodium hyaluronate, lidocaine hydrochloride, or to other amide-type local anesthetics.
- In pregnant or breast-feeding women
- In young people less than 18 years old

BELOTERO Balance Lidocaine is contra-indicated in patients who have a predisposition to keloid scarring.

Do not inject BELOTERO Balance Lidocaine into blood vessels.

Do not inject BELOTERO Balance Lidocaine to increase the volume of the breasts.

Do not inject BELOTERO Balance Lidocaine in the bones, tendons, ligaments, or muscles.

Do not inject BELOTERO Balance Lidocaine into the peri-orbital region (rings around the eyes, eyelids, crow's feet, wrinkles).

Do not use BELOTERO Balance Lidocaine to correct superficial wrinkles.

Do not inject BELOTERO Balance Lidocaine into areas presenting cutaneous problems of the inflammatory or infectious type, such as acne or herpes.

Do not use BELOTERO Balance Lidocaine in association with other aesthetic medicine techniques such as peeling, dermabrasion, or any type of laser treatment.

### Precautions for use

Before treatment, the patient must be informed about the device, its contra-indications and possible side effects.

Health care practitioners are encouraged to discuss all potential risks of soft tissue injection with their patients prior to treatment and ensure that patients are aware of signs and symptoms of potential complications.

Inadvertent intravascular injection, in particular in the glabellar area, could potentially lead to local vascular occlusion, ischemia and necrosis.

Patients suffering from rosacea should be informed that the treatment may increase the visibility of the rosacea.

In the absence of available clinical data on tolerance on the injection of BELOTERO Balance Lidocaine in patients with antecedents or with an active autoimmune disease or presenting history of severe multiple allergies or anaphylactic shock, the practitioner must decide whether to inject BELOTERO Balance Lidocaine on a case-by-case basis depending on the nature of the disease as well as the associated treatment. It is recommended to propose a prior double test on these patients and to not inject if the disease is active or evolving. It is also recommended to carefully monitor these patients after injection.

Patients with a history of streptococcal diseases or in patients pre-disposed to hypertrophic scars should also discuss about these conditions with the practitioner before proceeding with the aesthetic treatment.

Patients using anti-coagulant or anti-thrombolytic substances like aspirin or non-steroidal anti-inflammatory medications may have increased reactions of hematomas or bleeding at the injection site. Special care should apply in patients suffering from severe bleeding disorders.

In cases of patients suffering of epilepsy, impaired cardiac conditions, severely impaired hepatic function or severe renal dysfunction or porphyria, the practitioner must decide whether to inject BELOTERO Balance Lidocaine on a case-by-case basis depending of the nature of the disease as well as the associated treatment.

It has to be considered by practitioners and athletes that lidocaine may produce a positive results in anti-doping tests.

No clinical data is available on the injection of BELOTERO Balance Lidocaine into an area already treated with other aesthetic products or procedures. It is recommended not to inject into an area treated with a permanent or semi-permanent implant.

In order to minimize the risks of potential complications, this product should only be used by health care practitioners who have appropriate training, experience, and who are knowledgeable about the anatomy at and around the site of injection.

Check the integrity of the inner packaging prior to use, and the expiry date both for the syringe and the needle. Do not use these products if the expiry date has lapsed or if the packaging has been opened or damaged.

Needle has to be handled with care after removal of its protective shield.

If the needle becomes obstructed and the injection pressure becomes too high, stop the injection and change the needle or the cannula.

It should be noted that the presence of lidocaine may cause local redness or hypersensitivity.

Do not transfer BELOTERO Balance Lidocaine into another container and do not add other substances to the product.

Only the gel is sterile, but not the outside of the syringe. The outside walls of the syringe have not to be put in direct or indirect contact with the injection site or the needle.

BELOTERO Balance Lidocaine must not be used with an automated injection system.

Patient label is to be attached to patient records to ensure traceability of the product.

Do not recap the needle at the end of the injection session and discard it with the syringe and with the remaining product

after use. Disposal should be in accordance with the accepted medical practice and current applicable directives to ensure their correct elimination.

Do not re-sterilize and do not reuse, due to the associated risks including infection.

The patient should avoid applying makeup for at least 12 hours after treatment and avoid exposure to heat (saunas, Turkish baths, etc.), UV rays or massage and prolonged exposure to the sun for 2 weeks after the treatment. Sodium hyaluronate precipitates in the presence of quaternary ammonium salts (such as benzalkonium chloride). It is therefore recommended that BELOTERO Balance Lidocaine not come in contact with such substances.

There is no known interaction with other local or loco-regional anesthetics.

### Warnings

Introduction of product into the vasculature may lead to embolization, occlusion of the vessels, ischemia, or infarction

Rare but serious adverse events associated with the intravascular injection of soft tissue fillers in the face have been reported and include temporary or permanent vision impairment, blindness, cerebral ischemia or cerebral hemorrhage, leading to stroke, skin necrosis, and damage to underlying facial structures

Immediately stop the injection if a patient exhibits any of the following symptoms, including changes in vision, signs of a stroke, blanching of the skin, or unusual pain during or shortly after the procedure

Patients should receive prompt medical attention and possibly evaluation by an appropriate health care practitioner specialist should an intravascular injection occur.

### Side effects

Possible side effects exist and must be described to the patient by the practitioner before treatment. Very slight bleeding may occur during the injection although this disappears spontaneously as soon as the injection is finished. In occasional cases, one or more of the following may occur either immediately or as a delayed reaction (list not exhaustive):

- Reactions usually associated with injections like redness, erythema, edema, pain, burning sensation, sensitivity, sometimes accompanied by itching or stinging in the treated area. These reactions may last for a week,
- Hematomas in the treated area,
- Swelling in the treated area,
- Indurations or nodules in the treated area,
- Coloration or discoloration in the treated area,
- Injection site warmth, dryness, erosion,
- Injection site vesicles, pustules,
- Allergy to one of the product's components, especially sodium hyaluronate and lidocaine hydrochloride,
- Cases of necrosis, abscesses and granulomas after sodium hyaluronate injections have been reported in the literature. These rare potential risks must nevertheless be taken into account.

It is recommended that the patient report any undesirable effects lasting more than a week to his/her practitioner. The practitioner may then prescribe the patient appropriate treatment.

### Incident Reporting

Any incident that directly or indirectly led, might have led or might lead to any of the following: the death of a patient, user or other person, the temporary or permanent serious deterioration of a patient's, user's or other person's state of health, or a serious public health threat; and that has occurred in relation to the device should be reported under [AxCA-adverse.events@merz.com](mailto:AxCA-adverse.events@merz.com) or 1-877-336-4008 (Customer Service).

In case of serious incident, please contact in addition [Ax-Safety@merz.de](mailto:Ax-Safety@merz.de)

For Quality complaints please contact [AxCA-QA@merz.com](mailto:AxCA-QA@merz.com) or 1-877-336-4008 (Customer Service).

### Assembly of needle to syringe

For optimal use of BELOTERO Balance Lidocaine, it is important that the needle be properly connected to the syringe.

Improper assembly may result in separation of the needle and syringe during injection and/or product leakage at Luer-lock level. See diagrams 1, 2, 3 and 4.

1. **Firmly hold** the glass cylinder of the **syringe and the Luer-lock adaptor** between the thumb and the index fingers.
2. With the other hand, grasp the protective cap and unscrew it.
3. **Push & Twist** the needle on the syringe **until a resistance is felt**. Do not overtight. Overtight of the needle may lead to the Luer-lock moving and dislodging from the syringe.
4. Keep holding the Luer-lock and remove the sheath from the needle.

### Storage

Store between 2°C and 25°C. Protect from light and freezing. Avoid mechanical shocks.

### References

Updated documentation is available from ANTEIS SA, Switzerland.

Caution

Temperature limit of storage: 2 °C – 25 °C

Consult instruction for use  
(electronic instruction for use)

Date of manufacture

Do not use if package is damaged

Use-by date

Single use product. Do not re-use

Batch code

Open the blister by pulling the Tyvek lid following the arrow

Manufacturer

Sterile. Sterilized by moist heat. Only the gel is sterile, but not the outside of the syringe.







Sterile. Sterilized by irradiation. Only the needle itself is sterile, but not the outside of the needle packaging.

### Manufactured by:

ANTEIS SA  
18 Chemin des Aulx  
1228, Plan-les-Ouates  
Geneva/Switzerland  
reception@merz.com  
www.anteis.com

In case, you wish to have a paper instruction for use, please contact us  
phone 1-877-336-4008  
or email customerservicecan@merz.com

Date of revision of the instructions for use: 2026-01-01

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