Backstop in the right position during injection.

en INSTRUCTIONS FOR USE FOR BELOTERO® INTENSE LIDOCAINE

Description

BELOTERO Intense Lidocaine is a sterile, non-pyrogenic, viscoelastic, colourless, transparent gel of cross-linked sodium hyaluronate of non-animal origin, in a physiological phosphate buffer. BELOTERO Intense Lidocaine contains 0.3% of lidocaine hydrochloride.

Presentation

BELOTERO Intense Lidocaine is presented in a 1 ml sterile glass syringe for single use, sterilized in autoclave by moist heat. Each box contains two sterile CE-marked needles of 27G½ (0.4x13mm) and two traceability labels per syringe as well as instructions for use.

To ensure optimal use of BELOTERO Intense Lidocaine, it is recommended to use the needles provided in the box. The syringe is equipped with a Luer-lock and needle holding device that provides better handling and avoid disengagement of the needle and leakage at Luer-lock level during injection.

Composition

Cross-linked sodium hyaluronate: 25.5mg/ml

Lidocaine hydrochloride: 3.0 mg/ml

Phosphate buffer solution q.s.p.: 1.0 ml

The gel and all packaging items in direct contact with the gel are latex-free, bisphenol-A-free and phthalates-free.

The volume of the gel in each syringe is stated on the external box.

Indications

BELOTERO Intense Lidocaine is an injectable implant indicated to restore or significantly increase the volume of skin tissue which enables:

- Correction of facial depressions, either due to injury or age-related (scarring, moderate or deep wrinkles, furrows),
- Remodelling the contours of the face by injection into facial areas such as the cheeks, chin, etc.
- Increasing the volume of the lips (the border of the red part of the lips, or increasing the roundness)

The presence of lidocaine reduces local pain associated to the injection of the gel.

Posology and administration method

This device is designed to be injected into the deep dermis by a legally approved practitioner. For successful treatment it is essential that the practitioner has received prior training in the injection technique for soft tissue augmentation. The treatment must be carried out under the appropriate aseptic conditions. This product has been designed for use in physician's office. BELOTERO Intense Lidocaine must be injected into a healthy, non-inflamed, disinfected skin. To ensure optimal use of BELOTERO Intense Lidocaine it is recommended to assemble the needle according to the diagrams below. Start the injection pointing the needle bevel downwards.

The use of the included 27G½ needles is required, because a lower diameter would require greater force when injecting the implant.

BELOTERO Intense Lidocaine may be injected into the deep dermis using serial punctual injections, the linear retracing technique, anterotracing (push ahead) technique, cross-hatching technique (blanket) and the fan distribution technique.

BELOTERO Intense Lidocaine may be used in combination with BELOTERO Soft and BELOTERO Balance respecting the correct depth of injection and indication for each product.

If the needle becomes obstructed and the injection pressure becomes too high, stop the injection and change the needle.

Inject BELOTERO Intense Lidocaine slowly and apply the least amount of pressure necessary. The quantity to inject depends on the correction to be achieved. Over-correction is not necessary. After the injection, the practitioner can massage gently in order to distribute the product uniformly.

Contra-indications

BELOTERO Intense Lidocaine is contra-indicated in case of:

- known hypersensitivity to one of the product's components, especially sodium hyaluronate, lidocaine hydrochloride or to other amide-type local anesthetics.
- In pregnant or breast-feeding women
- In young people less than 18 years old

BELOTERO Intense Lidocaine is contra-indicated in patients who have a predisposition to keloid scarring.

Do not inject BELOTERO Intense Lidocaine into blood vessels.

Do not inject BELOTERO Intense Lidocaine to increase the volume of the breasts.

Do not inject BELOTERO Intense Lidocaine in the bones, tendons, ligaments, or muscles.

Do not inject BELOTERO Intense Lidocaine into the peri-orbital region (rings around the eyes, eyelids, crow's feet, wrinkles).

Do not use BELOTERO Intense Lidocaine to correct superficial wrinkles.

Do not inject BELOTERO Intense Lidocaine into areas presenting cutaneous problems of the inflammatory or infectious type, such as acne or herpes.

Do not use BELOTERO Intense Lidocaine in association with other aesthetic medicine techniques such as peeling, dermabrasion, or any type of laser treatment.

Precautions for use

Before treatment, the patient must be informed about the device, its contra-indications and possible side effects.

Health care practitioners are encouraged to discuss all potential risks of soft tissue injection with their patients prior to treatment and ensure that patients are aware of signs and symptoms of potential complications.

Inadvertent intravascular injection, in particular in the glabellar area, could potentially lead to local vascular occlusion, ischemia and necrosis.

Patients suffering from rosacea should be informed that the treatment may increase the visibility of the rosacea.

In the absence of available clinical data on tolerance on the injection of BELOTERO Intense Lidocaine in patients with antecedents or with an active autoimmune disease or presenting history of severe multiple allergies or anaphylactic shock, the practitioner must decide whether to inject BELOTERO Intense Lidocaine on a case-by-case basis depending on the nature of the disease as well as the associated treatment. It is recommended to propose a prior double test on these patients and to not inject if the disease is active or evolving. It is also recommended to carefully monitor these patients after injection.

Patients with a history of streptococcal diseases or in patients pre-disposed to hypertrophic scars should also discuss about these conditions with the practitioner before proceeding with the aesthetic treatment.

Patients using anti-coagulant or anti-thrombolytic substances like aspirin or non-steroidal anti-inflammatory medications may have increased reactions of hematomas or bleeding at the injection site. Special care should apply in patients suffering from severe bleeding disorders.

In cases of patients suffering of epilepsy, impaired cardiac conditions, severely impaired hepatic function or severe renal dysfunction or porphyria, the practitioner must decide whether to inject BELOTERO Intense Lidocaine on a case-by-case basis depending of the nature of the disease as well as the associated treatment.

It has to be considered by practitioners and athletes that lidocaine may produce a positive results in anti-doping tests.

No clinical data is available on the injection of BELOTERO Intense Lidocaine into an area already treated with other aesthetic products or procedures. It is recommended not to inject into an area treated with a permanent or semi-permanent implant.

In order to minimize the risks of potential complications, this product should only be used by health care practitioners who have appropriate training, experience, and who are knowledgeable about the anatomy at and around the site of injection.

Check the integrity of the inner packaging prior to use, and the expiry date both for the syringe and the needle. Do not use these products if the expiry date has lapsed or if the packaging has been opened or damaged.

Needle has to be handled with care after removal of its protective shield.

If the needle becomes obstructed and the injection pressure becomes too high, stop the injection and change the needle or the cannula.

It should be noted that the presence of lidocaine may cause local redness or hypersensitivity.

Do not transfer BELOTERO Intense Lidocaine into another container and do not add other substances to the product.

Only the gel is sterile, but not the outside of the syringe. The outside walls of the syringe have not to be put in direct or indi-

rect contact with the injection site or the needle.

BELOTERO Intense Lidocaine must not be used with an automated injection system.

Patient label is to be attached to patient records to ensure traceability of the product.

Do not recap the needle at the end of the injection session and discard it with the syringe and with the remaining product after use. Disposal should be in accordance with the accepted medical practice and current applicable directives to ensure their correct elimination.

Do not re-sterilize and do not reuse, due to the associated risks including infection.

The patient should avoid applying makeup for at least 12 hours after treatment and avoid exposure to heat (saunas, Turkish baths, etc.), UV rays or massage and prolonged exposure to the sun for 2 weeks after the treatment.

Sodium hyaluronate precipitates in the presence of quaternary ammonium salts (such as benzalkonium chloride). It is therefore recommended that BELOTERO Intense Lidocaine not come in contact with such substances.

There is no known interaction with other local or loco-regional anesthetics.

Warnings

Introduction of product into the vasculature may lead to embolization, occlusion of the vessels, ischemia, or infarction.

Rare but serious adverse events associated with the intravascular injection of soft tissue fillers in the face have been reported and include temporary or permanent vision impairment, blindness, cerebral ischemia or cerebral hemorrhage, leading to stroke, skin necrosis, and damage to underlying facial structures.

Immediately stop the injection if a patient exhibits any of the following symptoms, including changes in vision, signs of a stroke, blanching of the skin, or unusual pain during or shortly after the procedure.

Patients should receive prompt medical attention and possibly evaluation by an appropriate health care practitioner specialist should an intravascular injection occur.

Side effects

Possible side effects exist and must be described to the patient by the practitioner before treatment. Very slight bleeding may occur during the injection although this disappears spontaneously as soon as the injection is finished. In occasional cases, one or more of the following may occur either immediately or as a delayed reaction (list not exhaustive):

- Reactions usually associated with injections like redness, erythema, edema, pain, burning sensation, sensitivity, sometimes accompanied by itching or stinging in the treated area. These reactions may last for a week,
- Hematomas in the treated area,
- Swelling in the treated area,
- Indurations or nodules in the treated area,
- Coloration or discoloration in the treated area,
- Injection site warmth, dryness, erosion,
- Injection site vesicles, pustules,
- Allergy to one of the product's components, especially sodium hyaluronate and lidocaine hydrochloride,
- Cases of necrosis, abscesses and granulomas after sodium hyaluronate injections have been reported in the literature. These rare potential risks must nevertheless be taken into account.

It is recommended that the patient report any undesirable effects lasting more than a week to his/her practitioner. The practitioner may then prescribe the patient appropriate treatment.

Incident Reporting

Any incident that directly or indirectly led, might have led or might lead to any of the following: the death of a patient, user or other person, the temporary or permanent serious deterioration of a patient's, user's or other person's state of health, or a serious public health threat; and that has occurred in relation to the device should be reported under AxCA-adverse.events@merz.com or 1-877-336-4008 (Customer Service).

In case of serious incident, please contact in addition Ax-Safety@merz.de

For Quality complaints please contact AxCA-QA@merz.com or 1-877-336-4008 (Customer Service).

Assembly of needle to syringe

For optimal use of BELOTERO Intense Lidocaine, it is important that the needle be properly connected to the syringe.

Improper assembly may result in separation of the needle and syringe during injection and/or product leakage at Luer-lock level. See diagrams 1, 2, 3 and 4.

1. **Firmly hold** the glass cylinder of the **syringe and the Luer-lock adaptor** between the thumb and the index fingers.
2. With the other hand, grasp the protective cap and unscrew it.
3. **Push & Twist** the needle on the syringe **until a resistance is felt**. Do not overtight. Overtight of the needle may lead to the Luer-lock moving and dislodging from the syringe.
4. Keep holding the Luer-lock and remove the sheath from the needle.

Storage

Store between 2 °C and 25 °C. Protect from light and freezing. Avoid mechanical shocks.

References

Updated documentation is available from ANTEIS SA, Switzerland.

Caution

Temperature limit of storage: 2 °C – 25 °C

Consult instruction for use
(electronic instruction for use)

Date of manufacture

Do not use if package is damaged

Use-by date

Single use product. Do not re-use

Batch code

Open the blister by pulling the Tyvek lid following the arrow

Manufacturer

Sterile. Sterilized by moist heat. Only the gel is sterile, but not the outside of the syringe.







Sterile. Sterilized by irradiation. Only the needle itself is sterile, but not the outside of the needle packaging.

Manufactured by:

ANTEIS SA
18 Chemin des Aulx
1228, Plan-les-Ouates
Geneva/Switzerland
reception@merz.com
www.anteis.com

In case, you wish to have a paper instruction for use, please contact us
phone 1-877-336-4008
or email customerservicecan@merz.com

Date of revision of the instructions for use: 2026-01-01

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