Backstop in the right position during injection.

en INSTRUCTIONS FOR USE FOR BELOTERO® BALANCE LIDOCAINE

Description

BELOTERO Lips Shape is a sterile, non-pyrogenic, viscoelastic, colourless, transparent cross-linked sodium hyaluronate gel of non-animal origin in a physiological phosphate buffer. BELOTERO Lips Shape contains 0.3% of lidocaine hydrochloride.

Presentation

BELOTERO Lips Shape is presented in a single use pre-filled sterile glass syringe sterilized by moist heat. Each box contains one instruction leaflet, one syringe, two traceability labels and two sterile CE-marked needles for single use only. The dimensions of needles are stated on the external box.

Composition

Cross-linked sodium hyaluronate: 25.5 mg/ml

Lidocaine hydrochloride: 3.0 mg/ml

Phosphate buffer pH 7 q.s.: gel volume

The volume of the gel in each syringe is as stated on the external box.

Indications

BELOTERO Lips Shape is an injectable implant indicated to restore or significantly increase the volume of the skin tissue which enables increasing the volume of the lips (the border of the red part of the lips, or increasing the roundness) The presence of lidocaine reduces local pain associated to the injection of the gel.

Posology and administration method

BELOTERO Lips Shape is designed to be injected into the deep dermis and the mucous membrane of the lips by a legally approved practitioner.

BELOTERO Lips Shape can be used for all skin types.

For successful treatment it is essential that the practitioner has received a specific training on the injection techniques for soft tissue augmentation. A good knowledge of the anatomy and physiology of the site to be treated is required.

The treatment must be carried out under appropriate aseptic conditions.

BELOTERO Lips Shape must be injected into a healthy, non-inflamed and previously rigorously disinfected skin.

It is recommended to use one of the supplied needles.

To ensure optimal use of BELOTERO Lips Shape, it is recommended to assemble the needle according to the diagrams below.

The use of the enclosed 27G½ needle is recommended, as a smaller needle diameter would require a greater force to inject the implant.

BELOTERO Lips Shape can be used in combination with other BELOTERO products during the same session. Instructions for use of each product should be followed.

If the needle becomes obstructed and the injection pressure becomes too high, stop the injection and change the needle.

Inject BELOTERO Lips Shape slowly and apply the least amount of pressure necessary. The quantity of gel to be injected depends on the area to be treated and the correction to be achieved. Do not over-correct.

Gently massage the treated area after the injection to distribute the product uniformly.

Contra-indications

BELOTERO Lips Shape is contra-indicated:

- In case of known hypersensitivity to one of the product's components, especially to sodium hyaluronate, lidocaine hydrochloride and other amide-type local anesthetics,
- In pregnant or breast-feeding women,
- In young patients under 18 years old,
- In patients presenting a general infection,
- In patients who have a predisposition to keloid scarring.

Do not inject BELOTERO Lips Shape into blood vessels.

Do not inject BELOTERO Lips Shape into areas presenting cutaneous problems of an inflammatory or infectious type (acne, herpes...).

Do not inject BELOTERO Lips Shape into an area previously treated with a permanent dermal filler.

Do not inject BELOTERO Lips Shape in the glabellar region.

Do not inject BELOTERO Lips Shape into the periorbital region (rings around the eyes, eyelids, crow's feet, wrinkles).

Do not use BELOTERO Lips Shape to correct superficial wrinkles.

Precautions for use

Before treatment the patient must be informed about the device, its contra-indications and possible side effects.

Health care practitioners are encouraged to discuss all potential risks of soft tissue injection with their patients prior to treatment and ensure that patients are aware of signs and symptoms of potential complications

In the absence of available clinical data on tolerance and efficacy of the injection of BELOTERO Lips Shape in patients with antecedents or with an active autoimmune disease or in patients presenting a history of severe multiple allergies or anaphylactic shock, the practitioner must decide whether to inject BELOTERO Lips Shape on a case-by-case basis depending on the nature of the disease as well as the associated treatment. It is recommended to propose a prior double test to these patients and to not inject if the disease is evolving. It is also recommended to carefully monitor these patients after injection.

It is recommended not to inject BELOTERO Lips Shape in patients with a history of streptococcal diseases or in patients pre-disposed to hypertrophic scars or cheloids.

BELOTERO Lips Shape must not be used in association with other aesthetic medicine techniques such as peeling, dermabrasion, or any type of laser treatment before complete healing of the last treatment. In any case, even if the healing occurs earlier, BELOTERO Lips Shape must not be used earlier than 2 weeks after the last treatment. No clinical data is available on the combined use of BELOTERO Lips Shape with the above-mentioned treatments.

No clinical data is available on the injection of BELOTERO Lips Shape into an area already treated with other aesthetic products or procedures.

In order to minimize the risks of potential complications, this product should only be used by health care practitioners who have appropriate training, experience, and who are knowledgeable about the anatomy at and around the site of injection.

Check the integrity of the inner packaging and the expiry date for both the syringe and the needle prior to use. Do not use these products if the expiry date has lapsed or if the packaging has been opened or damaged.

Patients using antithrombotic substances such as aspirin or non-steroidal anti-inflammatory medications may have increased reactions of hematomas, nodules or bleeding at the injection site.

In cases of patients suffering from epilepsy, impaired cardiac conditions, severely impaired hepatic function or severe renal dysfunction or porphyria, the practitioner must decide whether to inject BELOTERO Lips Shape on a case-by-case basis depending of the nature of the disease as well as the associated treatment.

Practitioners and athletes should consider that lidocaine may produce positive results to anti-doping tests.

It should be noted that the presence of lidocaine may cause local redness or hypersensitivity.

For normal healthy adults, it is recommended that the maximum total dose of lidocaine HCl (without epinephrine) does not exceed 300 mg (or 4.5 mg/kg) per session. Overdosage of lidocaine HCl usually results in sign of the central nervous system or cardiovascular toxicity.

When using concurrently (topical administration...), the total administered dose of lidocaine should be considered. The concomitant use of other local anesthetic agents or agents structurally related to other amide-type local anesthetics should also be considered since the systemic toxic effects may be additive.

Care should be taken for patients with congenital methemoglobinemia, with glucose-6-phosphate dehydrogenase deficiencies and patients who are receiving concomitant treatment with methemoglobin-inducing agents.

Do not transfer BELOTERO Lips Shape into another container and do not add other substances to the product.

Only the gel is sterile, but not the outside of the syringe.

BELOTERO Lips Shape must not be used with an automated injection system not recommended by MERZ/ANTEIS. If an auto-

matic system is used, it is recommended that the practitioner has previously read its instructions for use and is trained on the use of the system.

Discard the syringe and the remaining product after use.

Do not re-sterilize and do not reuse due to the associated risks including infection.

The patient must avoid applying makeup for at least 12 hours after treatment as well as avoid saunas, Turkish baths and prolonged exposure to the sun or UV rays for 2 weeks after the treatment. Patients should also avoid putting pressure on and/or handling the treated area.

Sodium hyaluronate precipitates in the presence of quaternary ammonium salts (such as benzalkonium chloride). It is therefore recommended that BELOTERO Lips Shape does not come into contact with such substances.

There is no known interaction with other local or loco-regional anesthetics.

Warnings

Introduction of product into the vasculature may lead to embolization, occlusion of the vessels, ischemia, or infarction.

Rare but serious adverse events associated with the intravascular injection of soft tissue fillers in the face have been reported and include temporary or permanent vision impairment, blindness, cerebral ischemia or cerebral hemorrhage, leading to stroke, skin necrosis, and damage to underlying facial structures.

Immediately stop the injection if a patient exhibits any of the following symptoms, including changes in vision, signs of a stroke, blanching of the skin, or unusual pain during or shortly after the procedure

Patients should receive prompt medical attention and possibly evaluation by an appropriate health care practitioner specialist should an intravascular injection occur.

Side effects

Patients must be informed by the practitioner about possible side effects before treatment. A slight bleeding may occur during the injection, and it disappears spontaneously as soon as the injection is finished. In occasional cases one or more of the following may occur either immediately or as a delayed reaction (list not exhaustive):

- Reactions usually associated with injections such as redness, erythema, oedema or pain, sometimes accompanied by itching in the treated area. These reactions may last for a week.
- Hematomas in the treated area,
- Swelling in the treated area,
- Indurations or nodules in the treated area,
- Coloration or discoloration in the treated area,
- Injection site warmth, dryness, erosion,
- Injection site vesicles, pustules,
- Allergy to one of the product's components, especially to sodium hyaluronate and lidocaine hydrochloride.
- Cases of necrosis, abscesses and granulomas after sodium hyaluronate injections have been reported in the literature. These rare potential risks must nevertheless be considered.

Patients should be instructed to report any side effects which last for more than one week to his/her practitioner. The practitioner may then prescribe the patient appropriate treatment.

Incident Reporting

Any incident that directly or indirectly led, might have led or might lead to any of the following: the death of a patient, user or other person, the temporary or permanent serious deterioration of a patient's, user's or other person's state of health, or a serious public health threat; and that has occurred in relation to the device should be reported under AxCA-adverse.events@merz.com or 1-877-336-4008 (Customer Service).

In case of serious incident, please contact in addition Ax-Safety@merz.de

For Quality complaints please contact AxCA-QA@merz.com or 1-877-336-4008 (Customer Service).

Assembly of needle to syringe

For optimal use of Belotero Lips Shape, it is important that the needle is properly connected to the syringe. See diagrams 1, 2, 3 and 4.

1. **Firmly hold** the glass cylinder of the syringe and the Luer-lock adaptor between the thumb and the index fingers.
2. Grasp the protective cap with the other hand and unscrew it.
3. **Push & Twist** the needle on the syringe until a resistance is felt. Do not over-tight. Over-tightening of the needle may lead to the Luer-lock moving and dislodging from the syringe.
4. Keep holding the Luer-lock and remove the sheath from the needle.

Storage

Store between 2°C and 25°C. Protect from light and freezing. Avoid mechanical shocks.

References

Updated documentation may be available from ANTEIS SA, Switzerland.

Caution

Temperature limit of storage: 2 °C – 25 °C

Consult instruction for use
(electronic instruction for use)

Date of manufacture

Do not use if package is damaged

Use-by date

Single use product. Do not re-use

Batch code

Open the blister by pulling the Tyvek lid following the arrow

Manufacturer

Sterile. Sterilized by moist heat. Only the gel is sterile, but not the outside of the syringe.







Sterile. Sterilized by irradiation. Only the needle itself is sterile, but not the outside of the needle packaging.

Manufactured by:

ANTEIS SA
18 Chemin des Aulx
1228, Plan-les-Ouates
Geneva/Switzerland
reception@merz.com
www.anteis.com

In case, you wish to have a paper instruction for use, please contact us
phone 1-877-336-4008
or email customerservicecan@merz.com

Date of revision of the instructions for use: 2026-01-01

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