Backstop in the right position during injection.

en INSTRUCTIONS FOR USE FOR BELOTERO® VOLUME LIDOCAINE

Description

BELOTERO Volume Lidocaine is a sterile, non-pyrogenic, viscoelastic, colourless, transparent cross-linked sodium hyaluronate gel of non-animal origin, in a physiological phosphate buffer.

BELOTERO Volume Lidocaine contains 0.3% of lidocaine hydrochloride.

Presentation

BELOTERO Volume Lidocaine is presented in a single use pre-filled glass syringe, sterilized by moist heat.

Each box contains 2 syringes of BELOTERO Volume Lidocaine, 2 sterile 30G½" needles and 2 sterile 27G½" needles, 2 sterile cannulas 27G/37 mm, 1 instruction leaflet and 4 traceability labels.

The sterile CE-marked needles are for single use only.

Composition

Cross-linked sodium hyaluronate 26 mg/ml

Lidocaine hydrochloride 3,0 mg/ml

Phosphate buffer pH 7 q.s. as stated on the syringe label and the carton

The volume of the gel in each syringe is as stated on the syringe label and the carton.

Indications

BELOTERO Volume Lidocaine is an injectable resorbable implant indicated to restore facial volumes, as for example to enhance the cheeks or the chin.

The presence of lidocaine aims to reduce local pain associated with the injection of the gel and to improve patient comfort.

Posology and administration method

BELOTERO Volume Lidocaine is designed to be injected into the deep dermis, subcutaneously or in the upper periosteal, by a legally approved practitioner.

For successful treatment, it is essential that the practitioner has received specific training in injection techniques for restoring facial volume.

The treatment must be carried out under the appropriate aseptic conditions.

BELOTERO Volume Lidocaine must be injected into a healthy, non-inflamed and previously rigorously disinfected skin. To ensure optimal use of BELOTERO Volume Lidocaine, it is recommended to use the needles provided in the box. It is nonetheless possible to use other sterile needles or cannulae that can be adapted to the syringe's luer-lock connection depending on the injection technique chosen or the practitioner's preferences. In this case, it is possible to use a needle or a cannula recommended by Merz whose diameter is in the range 18G – 30G.

To ensure the optimal use of BELOTERO Volume Lidocaine, it is recommended to assemble the needle according to the diagrams below.

If the needle becomes obstructed and the injection pressure becomes too high, stop the injection and change the needle or the cannula.

Inject BELOTERO Volume Lidocaine slowly and apply the least amount of pressure necessary.

The quantity of gel to be injected depends on the area to be treated. It is recommended not to inject more than 3 ml per treatment site at each session. Do not over-correct.

After the injection, gently massage the treated area to distribute the product uniformly.

Contra-indications

BELOTERO Volume Lidocaine is contra-indicated:

- In case of known hypersensitivity to one of the product's components, especially sodium hyaluronate, lidocaine hydrochloride or to other amide-type local anesthetics,
- In pregnant and breast-feeding women,
- In young people less than 18 years old,
- In patients presenting a general infection,
- Predisposition to keloid scarring.

Do not inject BELOTERO Volume Lidocaine into blood vessels.

Do not inject BELOTERO Volume Lidocaine into the glabellar region.

Do not inject BELOTERO Volume Lidocaine into the tear-trough region

Do not inject BELOTERO Volume Lidocaine into the lips.

Do not inject BELOTERO Volume Lidocaine into areas presenting cutaneous problems of the inflammatory or infectious type (acne, herpes...).

Do not inject BELOTERO Volume Lidocaine for correction of superficial wrinkles and fine lines (injection of the product into the superficial dermis)

Do not inject BELOTERO Volume Lidocaine into an area previously treated with a permanent dermal filler

Do not inject BELOTERO Volume to increase the volume of the breasts.

Do not inject BELOTERO Volume in the bones, tendons and ligaments, or muscles.

Do not inject BELOTERO Volume into the periorbital region (rings around the eyes, eyelids, crow's feet, wrinkles).

BELOTERO Volume must not be used in association with other aesthetic medicine techniques such as peeling dermabrasion or any types of laser treatment.

Precautions for use

Before treatment, the patient must be informed about the device, its contra-indications and possible side effects.

Health care practitioners are encouraged to discuss all potential risks of soft tissue injection with their patients prior to treatment and ensure that patients are aware of signs and symptoms of potential complications.

In the absence of available clinical data on tolerance and efficacy on the injection of BELOTERO Volume Lidocaine in patients with antecedents or an active auto-immune disease, the practitioner must decide whether to inject BELOTERO Volume Lidocaine on a case-by-case basis depending on the nature of the disease as well as the associated treatment. It is recommended to propose a prior double test on these patients and to not inject if the disease is evolving. It is also recommended to carefully monitor these patients after injection.

BELOTERO Volume Lidocaine can be used to correct the signs of facial fat loss (lipoatrophy) in people with human immunodeficiency virus.

It is recommended not to inject BELOTERO Volume Lidocaine in patients with a history of streptococcal diseases or in patients pre-disposed to hypertrophic scars or cheloids.

It is recommended not to inject BELOTERO Volume Lidocaine into the periorbital region.

BELOTERO Volume Lidocaine must not be used in association with other aesthetic medicine techniques such as peeling, dermabrasion, or any type of laser treatment, and at the earliest after complete healing of the last treatment. In all cases, even if the healing occurs earlier, BELOTERO Volume Lidocaine must not be used earlier than 2 weeks after the last treatment. No clinical experiments are available on the combined use of BELOTERO Volume Lidocaine with the above-mentioned treatments.

No clinical data is available on the injection of BELOTERO Volume Lidocaine into an area already treated with other aesthetic products or procedures.

In order to minimize the risks of potential complications, this product should only be used by health care practitioners who have appropriate training, experience, and who are knowledgeable about the anatomy at and around the site of injection.

Check the integrity of the inner packaging prior to use, and the expiry date both for the syringe and the needle. Do not use these products if the expiry date has lapsed or if the packaging has been opened or damaged.

Patients using anti-coagulant substances like aspirin or non-steroidal anti-inflammatory medications may have increased reactions of haematomas or bleeding at the injection site.

In cases of patients suffering of epilepsy, impaired cardiac conditions, severely impaired hepatic function or severe renal dysfunction or porphyria, the practitioner must decide whether to inject BELOTERO Volume Lidocaine on a case-by-case basis depending of the nature of the disease as well as the associated treatment.

It has to be considered by practitioners and athletes that lidocaine may produce a positive results in anti-doping tests.

It should be noted that the presence of lidocaine may cause local redness or hypersensitivity.

For normal healthy adults, it is recommended that the maximum total dose of lidocaine HCl (without epinephrine) does not exceed 300 mg (or 4.5 mg/kg) per session. Overdosage of lidocaine HCl usually results in sign of the central nervous system or cardiovascular toxicity.

When using concurrently (topical administration...), the total administered dose of lidocaine should be considered. The concomitant use of other local anesthetic agents or agents structurally related to amide-type local anesthetics should also be considered since the systemic toxic effects may be additive.

Care should be taken for patients with congenital methemoglobinemia, with glucose-6-phosphate dehydrogenase deficiencies and patients who are receiving concomitant treatment with methemoglobin-inducing agents.

Do not transfer BELOTERO Volume Lidocaine into another container and do not add other substances to the product. Only the gel is sterile, but not the outside of the syringe.

BELOTERO Volume Lidocaine must not be used with an automated injection system.

Discard the syringe and the remaining product after use.

Do not re-sterilize and do not reuse, due to the associated risks including infection.

The patient must avoid applying makeup for at least 12 hours after treatment and avoid saunas and Turkish baths and prolonged exposure to the sun or UV rays for 2 weeks after the treatment. Patients should also avoid putting pressure on and/or handling the area of the face which has been treated.

Sodium hyaluronate precipitates in the presence of quaternary ammonium salts (such as benzalkonium chloride). It is therefore recommended that BELOTERO Volume Lidocaine not come in contact with such substances.

There is no known interaction with other local or loco-regional anesthetics.

Warnings

Introduction of product into the vasculature may lead to embolization, occlusion of the vessels, ischemia, or infarction.

Rare but serious adverse events associated with the intravascular injection of soft tissue fillers in the face have been reported and include temporary or permanent vision impairment, blindness, cerebral ischemia or cerebral hemorrhage, leading to stroke, skin necrosis, and damage to underlying facial structures.

Immediately stop the injection if a patient exhibits any of the following symptoms, including changes in vision, signs of a stroke, blanching of the skin, or unusual pain during or shortly after the procedure.

Patients should receive prompt medical attention and possibly evaluation by an appropriate health care practitioner specialist should an intravascular injection occur.

Side effects

Possible side effects exist and must be described to the patient by the practitioner before treatment. Very slight bleeding may occur during the injection although this disappears spontaneously as soon as the injection is finished. In occasional cases, one or more of the following may occur either immediately or as a delayed reaction (list not exhaustive):

- Reactions usually associated with injections like redness, erythema, oedema or pain sometimes accompanied by itching in the treated area. These reactions usually last for a week.
- Haematomas in the treated area.
- Swelling in the treated area.
- Indurations or nodules in the treated area.
- Coloration or discoloration in the treated area.
- Injection site warmth, dryness, erosion,
- Injection site vesicles, pustules,
- Allergy to one of the product's components, especially sodium hyaluronate and lidocaine hydrochloride.
- Cases of necrosis, abscesses and granulomas after sodium hyaluronate injections have been reported in the literature. These rare potential risks must nevertheless be taken into account.
- It is recommended that the patient report any undesirable effects lasting more than a week to his/her practitioner. The practitioner may then prescribe the patient appropriate treatment.

Incident Reporting

Any incident that directly or indirectly led, might have led or might lead to any of the following: the death of a patient, user or other person, the temporary or permanent serious deterioration of a patient's, user's or other person's state of health, or a serious public health threat; and that has occurred in relation to the device should be reported under AxCA-adverse events@merz.com or 1-877-336-4008 (Customer Service).

In case of serious incident, please contact in addition Ax-Safety@merz.de

For Quality complaints please contact AxCA-QA@merz.com or 1-877-336-4008 (Customer Service).

Assembly of needle to syringe

For optimal use of BELOTERO Volume Lidocaine, it is important that the needle/cannula be properly connected to the syringe. Improper assembly may result in separation of the needle/cannula and syringe during injection and/or product leakage at Luer-lock level. See diagrams 1, 2, 3 and 4.

1. Firmly hold the glass cylinder of the syringe and the Luer-lock adaptor between the thumb and the index fingers.
2. With the other hand, grasp the protective cap and unscrew it.
3. Push & Twist the needle/cannula on the syringe until a resistance is felt. Do not overtight
Overtight of the needle/cannula may lead to the Luer-lock moving and dislodging from the syringe.
4. Keep holding the Luer-lock and remove the sheath from the needle/cannula.

Storage

Store between 2 °C and 25 °C. Protect from light and freezing. Avoid mechanical shocks.

References

Updated documentation may be available from ANTEIS SA, Switzerland.

Caution

Temperature limit of storage: 2 °C – 25 °C

Consult instruction for use
(electronic instruction for use)

Date of manufacture

Do not use if package is damaged

Use-by date

Single use product. Do not re-use

Batch code

Open the blister by pulling the Tyvek lid following the arrow

Manufacturer

Sterile. Sterilized by moist heat. Only the gel is sterile, but not the outside of the syringe.







Sterile. Sterilized by irradiation. Only the needle itself is sterile, but not the outside of the needle packaging.

Manufactured by:

ANTEIS SA
18 Chemin des Aulx
1228, Plan-les-Ouates
Geneva/Switzerland
reception@merz.com
www.anteis.com

In case, you wish to have a paper instruction for use, please contact us
phone 1-877-336-4008
or email customerservicecan@merz.com

Date of revision of the instructions for use: 2026-01-01

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