
 Keim GmbH Rheinstraße 1-5, 63225 Langen Standort: Otto-Hahn-Straße 41, 63456 Hanau, Tel. +49 6181-9886-10	Job No.: me502744		Created at: 14.05.2018		Operator: bka
	Size: 297 x 210 mm (b x h)		1. AC: 02.07.2018 to	5. AC:	
	MN 62245	MC 0000	2. AC:	6. AC:	
	AM 31210	MZ/AZ 0000	3. AC:	7. AC:	
	Typestyle: Avenir		4. AC:	8. AC:	
Belotero GA Soft Lidocaine Kanada		 Black Stanzfarbe			

Attention
Caution

Consulter la notice
Consult instructions for use

Ne pas utiliser si l'emballage est endommagé
Do not use if package is damaged

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GB INSTRUCTIONS FOR USE FOR BELOTERO® SOFT LIDOCAINE	
Description	BELOTERO Soft Lidocaine is a sterile, non-pyrogenic, viscoelastic, colourless, transparent gel of cross-linked sodium hyaluronate of non-animal origin, in a physiological phosphate buffer. BELOTERO Soft Lidocaine contains 0.3% of lidocaine hydrochloride.
Presentation	BELOTERO Soft Lidocaine is presented in a 1 ml sterile glass syringe for single use, sterilized by moist heat. Each box contains two sterile CE-marked needles of 30G½ (0.3x13mm) and two traceability labels per syringe as well as instructions for use. To ensure optimal use of BELOTERO Soft Lidocaine, it is recommended to use the needles provided in the box.
	The syringe is equipped with a Luer-lock and needle holding device that provides better handling and avoid disengagement of the needle and leakage at Luer-lock level during injection.
Composition	Cross-linked sodium hyaluronate: 20.0 mg/ml Lidocaine hydrochloride: 3.0 mg/ml Phosphate buffer solution q.s.p.: 1.0 ml
	The gel and all packaging items in direct contact with the gel are latex-free, bisphenol-A-free and phthalates-free.
Indications	BELOTERO Soft Lidocaine is indicated for use to increase the volume of skin tissue which enables:
	• Correcting shallow facial depressions, either of accidental origin, or age-related (scars, fine lines, crow's feet wrinkles),
	• Retouching a correction undertaken with BELOTERO Soft Lidocaine, immediately or a few weeks after the first injection.
	The presence of lidocaine reduces local pain associated to the injection of the gel.
Posology and administration method	This device is designed to be injected into the superficial to mid-dermis by a legally approved practitioner. For successful treatment it is essential that the practitioner has received prior training in the injection technique for increasing the volume of skin tissues. The treatment must be carried out under the appropriate aseptic conditions. This product has been designed for use in physician's office. BELOTERO Soft Lidocaine must be injected into a healthy, non-inflamed, disinfected skin. To ensure optimal use of BELOTERO Soft Lidocaine, it is recommended to assemble the needle according to the diagrams below. Start the injection pointing the needle bevel downwards.
	BELOTERO Soft Lidocaine may be injected into the superficial dermis using serial punctual injections, the linear retracing technique, anterotracing (push ahead) technique, cross-hatching technique (blanket), fan distribution technique, blanching technique, multi-puncture and micro-papular injection techniques. BELOTERO Soft Lidocaine may be used in combination with BELOTERO Balance and BELOTERO Intense respecting the correct depth of injection and indication for each product. It can also be combined with the BELOTERO Hydro for the dermis rehydration purposes.
	If the needle becomes obstructed and the injection pressure becomes too high, stop the injection and change the needle. Inject BELOTERO Soft Lidocaine slowly and apply the least amount of pressure necessary. Over-correction is not necessary. After the injection, the practitioner can massage gently in order to distribute the product uniformly.
Contra-indications	BELOTERO Soft Lidocaine is contra-indicated in case of:
	• known hypersensitivity to one of the product's components, especially sodium hyaluronate, lidocaine hydrochloride or to other amide-type local anesthetics.
	• In pregnant or breast-feeding women
	• In young people less than 18 years old

BELOTERO Soft Lidocaine is contra-indicated in patients who have a predisposition to keloid scarring.

Do not inject BELOTERO Soft Lidocaine into blood vessels.

Do not inject BELOTERO Soft Lidocaine to increase the volume of the breasts. Do not inject BELOTERO Soft Lidocaine in the bones, tendons, ligaments, or muscles.

Do not inject BELOTERO Soft Lidocaine into the peri-orbital region (rings around the eyes, eyelids, crow's feet, wrinkles).

Do not inject BELOTERO Soft Lidocaine into areas presenting cutaneous problems of the inflammatory or infectious type, such as acne or herpes.

Do not use BELOTERO Soft Lidocaine in association with other aesthetic medicine techniques such as peeling, dermabrasion, or any type of laser treatment.

Precautions for use
Before treatment, the patient must be informed about the device, its contra-indications and possible side effects.

Health care practitioners are encouraged to discuss all potential risks of soft tissue injection with their patients prior to treatment and ensure that patients are aware of signs and symptoms of potential complications.

Inadvertent intravascular injection, in particular in the glabellar area, could potentially lead to local vascular occlusion, ischemia and necrosis.

Patients suffering from rosacea should be informed that the treatment may increase the visibility of the rosacea.

In the absence of available clinical data on tolerance on the injection of BELOTERO Soft Lidocaine in patients with antecedents or with an active autoimmune disease or presenting history of severe multiple allergies or anaphylactic shock, the practitioner must decide whether to inject BELOTERO Soft Lidocaine on a case-by-case basis depending on the nature of the disease as well as the associated treatment. It is recommended to propose a prior double test on these patients and to not inject if the disease is active or evolving. It is also recommended to carefully monitor these patients after injection.

Patients with a history of streptococcal diseases or in patients pre-disposed to hypertrophic scars should also discuss about these conditions with the practitioner before proceeding with the aesthetic treatment.

Patients using anti-coagulant or anti-thrombotic substances like aspirin or non-steroidal anti-inflammatory medications may have increased reactions of hematomas or bleeding at the injection site. Special care should apply in patients suffering from severe bleeding disorders.

In cases of patients suffering of epilepsy, impaired cardiac conditions, severely impaired hepatic function or severe renal dysfunction or porphyria, the practitioner must decide whether to inject BELOTERO products on a case-by-case basis depending of the nature of the disease as well as the associated treatment.

It has to be considered by practitioners and athletes that lidocaine may produce a positive results in anti-doping tests.

No clinical data is available on the injection of BELOTERO Soft Lidocaine into an area already treated with other aesthetic products or procedures. It is recommended not to inject into an area treated with a permanent or semi-permanent implant.

In order to minimize the risks of potential complications, this product should only be used by health care practitioners who have appropriate training, experience, and who are knowledgeable about the anatomy at and around the site of injection.

Check the integrity of the inner packaging prior to use, and the expiry date both for the syringe and the needle. Do not use these products if the expiry date has lapsed or if the packaging has been opened or damaged.

Needle has to be handled with care after removal of its protective shield. If the needle becomes obstructed and the injection pressure becomes too high, stop the injection and change the needle or the cannula.

It should be noted that the presence of lidocaine may cause local redness or hypersensitivity.

Do not transfer BELOTERO Soft Lidocaine into another container and do not add other substances to the product.

Only the gel is sterile, but not the outside of the syringe. The outside walls of the syringe have not to be put in direct or indirect contact with the injection site or the needle.

BELOTERO Soft Lidocaine must not be used with an automated injection system.

Patient label is to be attached to patient records to ensure traceability of the product.

Do not recap the needle at the end of the injection session and discard it with the syringe and with the remaining product after use. Disposal should be in accordance with the accepted medical practice and current applicable directives to ensure their correct elimination.

Do not re-sterilize and do not reuse, due to the associated risks including infection. The patient should avoid applying makeup for at least 12 hours after treatment and avoid exposure to heat (saunas, Turkish baths, etc.), UV rays or massage and prolonged exposure to the sun for 2 weeks after the treatment.

Sodium hyaluronate precipitates in the presence of quaternary ammonium salts (such as benzalkonium chloride). It is therefore recommended that BELOTERO Soft Lidocaine not come in contact with such substances.

There is no known interaction with other local or loco-regional anesthetics.

Warnings
Introduction of the product into the vasculature may lead to embolization, occlusion of the vessels, ischemia, or infarction.

Rare but serious adverse events associated with the intravascular injection of soft tissue fillers in the face have been reported and include temporary or permanent vision impairment, blindness, cerebral ischemia or cerebral hemorrhage, leading to stroke, skin necrosis, and damage to underlying facial structures.

Immediately stop the injection if a patient exhibits any of the following symptoms, including changes in vision, signs of a stroke, blanching of the skin, or unusual pain during or shortly after the procedure.

Patients should receive prompt medical attention and possibly evaluation by an appropriate health care practitioner specialist should an intravascular injection occur.

Side effects
Possible side effects exist and must be described to the patient by the practitioner before treatment. Very slight bleeding may occur during the injection although this disappears spontaneously as soon as the injection is finished. In occasional cases, one or more of the following may occur either immediately or as a delayed reaction (list not exhaustive):

- Reactions usually associated with injections like redness, erythema, edema, pain, burning sensation, sensitivity, sometimes accompanied by itching or stinging in the treated area. These reactions may last for a week,
- Hematomas in the treated area,
- Swelling in the treated area,
- Indurations or nodules in the treated area,
- Coloration or discoloration in the treated area,
- Allergy to one of the product's components, especially sodium hyaluronate and lidocaine hydrochloride,
- Cases of necrosis, abscesses and granulomas after sodium hyaluronate injections have been reported in the literature. These rare potential risks must nevertheless be taken into account.

It is recommended that the patient report any undesirable effects lasting more than a week to his/her practitioner. The practitioner may then prescribe the patient appropriate treatment.

Assembly of needle to syringe
For optimal use of BELOTERO Soft Lidocaine, it is important that the needle be properly connected to the syringe. Improper assembly may result in separation of the needle and syringe during injection and/or product leakage at Luer-lock level. See diagrams 1, 2, 3 and 4.

1. **Firmly hold the glass cylinder of the syringe and the Luer-lock adaptor** between the thumb and the index fingers.

2. With the other hand, grasp the protective cap and unscrew it.

3. **Push & Twist** the needle on the syringe until a resistance is felt. Do not over-tight. Overtight of the needle may lead to the Luer-lock moving and dislodging from the syringe.

4. Keep holding the Luer-lock and remove the sheath from the needle.

Storage
Store between 2°C and 25°C. Protect from light and freezing. Avoid mechanical shocks.

References
Updated documentation is available from ANTEIS SA, Switzerland.

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